

BRIEF DESCRIPTION of SEVERAL VIEWS of the DRAWING

0001 Figure 1 shows a view seen from the top of the saliva ejector tip. Seen are six holes of one mm diameter each, placed circumferentially, 2mm from the edge.

0002 Figure 2 shows a view seen from the side of the saliva ejector tip. Seen are four holes of one mm diameter each. The top two holes are two of six holes evenly spaced and placed circumferentially 3mm from the end of the tip. The lower two holes are two of six holes, evenly spaced and placed circumferentially 5mm from the end of the tip.